

Hands-On – Primers design for PCR

Objective: Design PCR primers and check them for specificity

Exercise 1:

Sequence1

```
TCATAGCTCCGATCCGTTCCGTTCAATCAATGTACCTGGCGGCCATTTCGTTCAACGGACGTCCAACCGGG
AATATCACGACGAGTGGGTCTCGGTGCGCGTGGAAGGGAGATCGGCGTTCACCTCTGCGATGCTCGACC
GCACTCGAACCCCCGTCAAATGCTAGCCATAACCGTTTTGGCAATATCTCTCAACGTGGTTCCGCATGGC
GTGCCGAACTGTCCCAAAGGCGAATGGGACACCAAGCTCCCGCAAAGTTGGAGTCAAAGCACCGACCGAG
TGAGCGATGCCCATTCGCTCCAGATCGGGCCGACCGCGTATAAGCCCAGTTCAGCCACAAGGAGATCAG
TCTACCAACCTTAATGAAACGGCGCAACACGCCCATGTGGCTTGCTATCCCGACCGAGTCGTAAGCAAT
TGCGCGGCGTTCGGTCTCGCGCCGGCCCAACTGCGGCCACCCTCAAATGGTTTTCGCGTGGAACGCTCCT
GTGGGCCCATAGGATTTTCGAAAAAATCTGAGAGCTCCTGGTGGTCTGCACGTTCCAGCTGATTTTCCC
AGCATAGCTTCCACTGCACTTTTAAGGACAT
```

Procedure:

- Use the given sequence to design a couple of primers for PCR.
- Go to the primer3 program: <http://frodo.wi.mit.edu/>
- Paste your sequence in the window
- In the options: Pick the left and right primers
- And you have to choose some options:
 - The product size ranges: 100–300 (you can give different ranges)
 - Specify the number of primers to return (you can choose as many as you want)
 - Optimum of T_m : 60°C

- Optimum of GC content: 50%

Next step: Check the primers for their specificity

We will perform an *Insilico* PCR amplification

- Use the following link (<http://insilico.ehu.es/PCR/>)
- Paste the picked primers in the corresponding window
- Check if the primers can amplify all the species of the *sinorhizobium*
- Specify the maximum length of bands
- Click the amplify button
- Interpret the result
- Get the Fasta format of the amplicons
- Check if the sequences of the amplicons are the same

Restriction mapping

In order to check if some sequences can be identified based on their restriction profile; we are going to perform an *insilico* restriction digest of complete genome:

- Will try the *sinorhizobium* genome
- Specify the “Only restriction enzymes with known bases (no N,R,Y...)”
- Click the restriction enzyme digest of genome
- Interpret the result
- Click the link in the right of the table for the PFGE simulation

Exercise 2:

The accession number FJ236988.1 correspond to eukaryotic gene.

We would like to design primers to amplify only mRNA, but have to flank the gene intron.

- Why is it interesting to have primers flanking intron?
- Use the primer3 program